
# Skandika

## HEIMTRAINER ERGOMETER ELSKLING E 520

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
# Table of contents

Important safety notice .....	27
Assembly instructions.....	28
Warm-up and cool-down.....	37
Computer operation manual.....	38
Using a chest-belt.....	42
Using the apps „Fitness Data“ & „KinoMap“ .....	43
Exploded diagram .....	46
Parts list.....	48
Cleaning and care .....	48
Guarantee Conditions.....	49

Visit our website  
for further information  
[www.skandika.com](http://www.skandika.com)

Scan the QR code  
with your smartphone

## **WARNING**

To reduce risk of injury, read and understand this instruction manual before using the device!  
This machine is intended for home use only in accordance with the instructions provided in this manual. Read the instruction manual carefully before using this device and keep the instruction manual for future use.






### **IMPORTANT SAFETY NOTICE**

- To ensure the best safety of the device, regularly check it on damages and worn parts.
- If you pass on this device to another person or if you allow another person to use it, make sure that that person is familiar with the content and instructions in these instructions.
- Only one person should use the exerciser at a time. The maximum load of this device is 135 kg (298 lbs). Class: H - not suitable for therapeutically use!
- Before the first use and regularly make sure that all screws, bolts and other joints are properly tightened and firmly seated.
- Before you start your workout, remove all sharp-edged objects around the equipment.
- Only use the exercise for your work-out if it works flawlessly.
- Any broken, worn or defective part must immediately be replaced and/or the device must no longer be used until it has been properly maintained and repaired. DO NOT use attachments not recommended by the manufacturer.
- Keep children and pets away from the machine at all times. DO NOT leave children unattended in the same room with the machine.
- Parents and other supervisory persons should be aware of their responsibility, due to situations which may arise for which the device has not been designed and which may occur due to children's natural play instinct and interest in experimenting.
- If you do allow children to use this device, be sure to take into consideration and assess their mental and physical condition and development, and above all their temperament. Children should use the device only under adult supervision and be instructed on the correct and proper use of the device. The device is not a toy.
- Make sure there is sufficient free space around the appliance when you set it up.
- Please note that an improper and excessive work-out may be harmful to your health.
- Please note that levers and other adjustment mechanisms are not projecting into the area of movement during the workout.
- When setting up the device, please make sure that the device is standing in a stable way and that any possible unevenness of the floor is evened out.
- Always wear appropriate clothing and shoes which are suitable for your workout on the device. The clothes must be designed in a way so that they will not get caught in any part of the device during the work-out due to their form (for example, length). Be sure to wear appropriate shoes which are suitable for the work-out, firmly support the feet and which are provided with a non-slip sole.
- Be sure to consult a physician before you start any device program. He may give you proper hints and advice with respect to the individual intensity of stress for you, your workout and sensible eating habits.
- The general rule is that devices and training devices are no toys. Therefore, they must only be used by properly informed or instructed persons.

## IMPORTANT SAFETY NOTICE

- Be sure to set up the device in a dry and even place and always protect it from humidity. If you wish to protect the place particularly against pressure points, contamination, etc., it is recommended to put a suitable, non-slip mat under the device. DO NOT use the machine near liquids or outdoors.
- Stop your workout immediately in case of dizziness, nausea, chest pain or any other physical symptoms. In case of doubt, consult your physician immediately.
- Before beginning any exercise program, consult your physician. This is especially important for individuals with pre-existing health problems.
- Children, disabled and handicapped persons should use the device only under supervision and in presence of another person who may give support and useful instructions.
- Be sure that your body parts and those of other persons are never close to any moving parts of the device during its use.
- When adjusting the adjustable parts, make sure they are adjusted properly and note the marked, maximum adjusting position, for example of the saddle support, respectively. Pay attention to the warnings on the device.
- Do not work out immediately after meals!
- Before using the device, always do stretching exercises to properly warm up.

## Assembly instructions - Small parts bag:

No.	Description	Q'ty		
59	Pan bolt M8 x 55	4		
71	Pan bolt M8 x 16	6		
70	Flat washer	2		
51	Screw	2		
22	Release knob	1		
57	Flat washer	8		
16	T-Knob	1		
46	Bushing	1		

# Assembly instructions - Checklist:

No.	Description	Q'ty		
1/25L/25R /27/37	Main frame / Chain cover left & right / Crank	1 Set		
5/50	Handlebars / Hand pulse	1 Set		
3/4/19/ 20/22	Seat tube / Slide bracket / Knob	1 Set		
2/43/ 48/49	Front post / Cover / Sensor cable	1 Set		
9/48/69	Computer with screws	1 Set		
8/56/57	OPTIONAL: Seat with fixing material	1 Set		
13	Front clamp cover for front post	1		
14	Rear clamp cover for front post	1		
7/29	Rear stabilizer w/ adjustable endcaps	1 Set		
6/40	Front stabilizer w/ transportation wheels	1 Set		
23L/23R	Pedal left/right	1 Set		
22	Release knob	1		
73	Adaptor	1		
72	Bottle holder	1		


Before starting to assemble we recommend to look for a place with sufficient space. Use a floor protection. Remove all parts from the carton and spread them on the floor. Check according to the checklist, if all parts are complete. Perhaps some parts are already pre-assembled for your convenience.

We recommend to set-up the device with 2 persons.

### Step 1

Fix the front stabilizer (6) on front part of main frame (1) with 2 screws (59) and 2 flat washers (57) as shown.

Then fix the rear stabilizer (7) on rear part of main frame (1) with 2 screws (59) and 2 flat washers (57).


## Step 2

Fix the covers (29) and (40) on their respective stabilizer with the screws (52) as shown.

Mount the right pedal (23R) onto the right crank arm and tighten it by turning it clockwise. Afterwards fix the left pedal (23L) accordingly on the left crank arm (to be screwed in anti-clockwise direction).


The pedals are marked with "R" (right) and "L" (left).

Fix the pedal strap on the pedals and adjust the strap according to your feet/shoe size.

### Step 3


Slide the cover (20) from below through the seat tube (3) as shown and insert it into the main frame (1). Align the drillings and fix the seat tube with the knob (22) at a proper height. The seat height may be adjusted according to your individual needs later on.

## Step 4

Fix the optional seat (8) on the sliding bracket (4) with the knob (19) and the flat washer (57).


## Step 5

Remove the pre-assembled 6 screws (71) and flat washers (57) from the main frame (1).

Slide the cover (43) from below through the front post (2) and connect the upper sensor cable (48) to the lower sensor cable (49).


Insert the front post (2) into the main frame (1) and fix it with the screws (71) and flat washers (57). Mount the bottle holder (72) on front post (2) with 2 screws (52) .

## Step 6

Fix the handlebars (5) on the front post (2) with the help of the T-knob (16), the front clamp cover (14), the bushing (46) and the screw (52).


Afterwards fix the rear clamp cover (13) with each 2 screws (51) and (52).

## Step 7

Loosen the 4 screws (69) from the back of computer (9). Connect the upper sensor cable (48) and the pulse cables (50) to the computer (9). Fix the computer (9) to the bracket of front post (2) using the 4 screws (69).

CHECK IF ALL BOLTS AND NUTS ARE TIGHTENED BEFORE USING THE DEVICE!

### Mains connection

The computer gets its power supply by the included adaptor. The adaptor will be connected to a normal mains plug (European standard, please check the specifications resp. please use a fitting mains adaptor for your power system). If you have to buy a new adaptor it is important to pay attention that it is GS approved and corresponds to the product safety law.

Don't use another electricity supply because this could cause damages or impairments of functions.

## WARM-UP and COOL-DOWN


A successful exercise program consists of a warm-up, aerobic exercise and a cool-down. Exercise for at least two or preferably three times a week, resting for a day between workouts. After several months, you may increase the frequency to four or five times per week.

### WARM-UP

This stage helps to improve the blood circulation and prepares the muscles for your workout. It additionally helps to reduce the risk of injury or cramps. It is recommended to do some stretching exercises as shown here below. Hold each stretching position for approximately 30 seconds. Never force or jerk yourself into a stretching position – if you feel pain, STOP immediately. Warm-up exercises may also include brisk walking, jogging, jumping jacks, jump rope exercises or running in place.

### STRETCHING

Muscles can be stretched more easily when these are warm. This reduces the risk of injury. DO NOT BOUNCE.


Remember always to check with your physician before starting any exercise program.

### COOL-DOWN

This stage helps to soothe your muscles and your cardiovascular system after your workout.

At the end of your workout, reduce the speed and exercise for approximately 5 minutes at this lower speed level. Afterwards, repeat the warm-up exercises as described above.

# Computer operation manual

## The meaning of the buttons

**Mode** : Used for choosing functions resp. programs.

**Recovery** : Press this button to start the heart recovery measurement.

**ST / SP** : 1) Start or finish a workout.

2) Start or finish a body fat measurement.

3) Press & hold this button for approx. 3 seconds to reset all values to 0.

**Enter** : Confirm setting or selection.

**Down (Ab) / Up (Auf)** : Decrease / Increase the resistance during exercise or adjust values.

## The computer functions

### TIME

Computer shows actual exercise time. Preset a target time: You may preset a target time (max. 99:00) before starting a workout. The time will then count down and an acoustic signal will be emitted when reaching 0:00.

### RPM

The computer shows the current rounds/min.in the display.

### SPEED

The computer shows the current speed in the display.

### DIST (Distance)

The workout distance is displayed here. Preset a target distance: You may preset a target distance (max. 99.9 km) before starting a workout. The distance will then count down and an acoustic signal will be emitted when reaching 0.0.

**ODO**

Shows the accumulated distance of all workouts with this device.

**CAL (Calories)**

The currently burned calories are displayed here. The calorie consumption values registered and displayed with this computer merely serve as indicative values for a person of medium stature with medium resistance setting, and can deviate considerably from the actual medically precise calorie consumption. You may also set a target calorie consumption before training, max. value = 999 kcal.

**WATT**

Computer shows actual watt value.

**PULSE**

The computer is showing the user's heart rate in beats per minute. For a correct display please ensure that you fully grasp the hand sensors. It may last up to 2 minutes before the pulse frequency will be measured correctly. You may also use a compatible chest belt - please find more about this topic in chapter „Using a chest belt“. If you use a chest belt, the measured frequency will also be displayed here.

## The different programs / modes

**P1 - Manual mode**

In manual mode you can adjust the load manually. You may select the tension level by using „UP“ and „DOWN“. Press „ST/SP“ to start or finish a workout. You may also enter target values for time and other functions.

**P2 to P14: Mountain- and valley courses**

Within this mode you can exercise with a preset programme profile and the computer will automatically change the tension level. You may additionally adjust the tension level with „UP“ and „DOWN“. You may also input target values for different function values (e.g. time etc.)

**P15: Exercise with preset Watt value**

In this exercise mode, you can preset WATT and other values. After you have selected this mode, the device will prompt you to input values (e.g. TIME, WATT etc.). Press „UP“ / „DOWN“ and „ENTER“ to input all these values. Press „ST/SP“ to start the exercise.

The Watt load depends on the rounds per minute and the level of resistance. The device is able to calculate the Watt load and display the respective value. Your computer will control the tension automatically during exercise. You can not adjust the resistance manually as usual in this programme!

Please pay attention, that if a watt preset value is of medical necessity in your case, you need definitely talk to your doctor before using this device. Your doctor should accompany this exercise and should give you advice on the use and surveillance.

**P16: Body-fat measurement**

This mode will process a body-fat measurement. In the following we will present you this technology and for what it is good. However, we strongly recommend to consult your physician before starting this body fat measurement in order to clarify if such a measurement by electrical current is suitable for you.

How to reach your ideal body type and improve your health

Being fit has many advantages, from helping you to control your weight to giving you a better night sleep. Most important of all, there is impressive and mounting evidence that people who exercise regularly, and in the right way, are less prone to killer diseases such as heart attacks and strokes and live longer than people who do not exercise.

Different diseases like e.g. diabetes, heart diseases and certain forms of cancer can be controlled or even avoided by controlling body fat levels.

If you watch your body fat you can reduce or increase your weight according to measured values by taking in more or less calories (ask a specialist for more information). Usually the individual part of body fat will be given as a percentage of the total body weight.

Weight reduction should happen by reduction of body fat and not by reduction of muscle mass or water.

#### Bio-electrical Impedance Analysis – Approved by Henry C. Lukaski, Ph. D. Grand Forks Human Nutrition Research Center

The fat monitoring determines body fat using Bio-Electric Impedance Analysis (BIA), a process by which a small electrical current is sent through the body to measure the impedance. The current will flow through the hands. The more fat containing adipose tissue (fat) in the body, the harder the current must work to travel around. In this case, the resistance to the current generates a higher impedance value, which in turn reflects a higher value for body fat. The impedance value is then inserted into an equation (along with height, weight, sex, and physical stature) which calculates body density and percent body fat.

#### Additional hints for your health

The body fat analysis is suitable for adults and children of 7 years age and older. This measurement is in no way suitable for people with heart diseases or when taking cardiovascular medicine. Please also do not use this item if you are pregnant. Wrong values can be measured if you belong to these groups:

- People over 70 years of age
- Dialyse-Patients
- People, who suffer from Oedema
- Intensive sport people, Bodybuilders, Professional athletes
- If when relaxing your heart rate is below 60

In these cases the result of the measurement cannot be regarded as correct value, but only as a guide for the improvement between different measurements. Always consult your physician before starting a body fat measurement! In case of improper handling or current fluctuations measured values can be lost and/or be displayed wrongly.

#### **How to do the bodyfat measurement**

The correct posture for the body fat measurement:

Sit on the bike in an upright position and do not bend your elbows and knees too strongly. Do not move during the measurement and adjust your seating height in the way, that the arms are pointing straight forward in the height of your heart while your hands are grasping the handpulse-sensors.

This is to ensure, that the current can flow in an optimal way through your body. After choosing the correct programme (**P16**) adjust the following values with the „UP“, „DOWN“ and „ENTER“ buttons:

1. Sex („GENDER“; Male or Female)
2. Height
3. Weight
4. Age

Take the correct posture for the measurement (as described earlier) and press „ST/SP“. During measuring, hold both hands firmly on the handsensors. The measurement will take several seconds.

After the measurement has ended the computer will display your bodyfat percentage (Fat%) as well as the BMI (Body-Mass-Index), BMR (Basal Metabolism Ration) and your body type. There are 6 different body types (depending on the measure bodyfat percentage) as follows:

Type 1: Men from 7 - 9,9%; Women from 14 - 16,9%; Type 2: Men from 10 - 12,9%; Women from 17 - 19,9%

Type 3: Men from 13 - 16,9%; Women from 20 - 23,9%; Type 4: Men from 17 - 19,9%; Women from 24 - 27,9%

Type 5: Men from 20 - 24,9%; Women from 28 - 29,9%; Type 6: Men > 24,9%; Women > 29,9%

Press ST/SP to leave the bodyfat mode.

#### *Hints for the shown results:*

Almost not one measurement of body fat will show the same results like before. Please keep in mind that (similar to blood pressure measurement) the smallest movements or surrounding influences or also your daily constitution will have a strong influence upon the results. So, it is normal if you have different results for one person within several measurements. Therefore the shown results are always only indications to give you hints for your body constitution and will help you to find the most effective way of training your body. The displayed values should only be regarded as evidence of your physical condition in order to support you finding your individual best way of exercise.

#### **P17 to P20 - Target heart rate training**

These modes require the input of a target heart rate. The device will then try to control the resistance level automatically according to the desired heart rate level. The programme 17 allows the direct input of a target heart rate - the programmes 18 to 21 will calculate the maximum pulse frequency by the input of your age and a percentage. The formula is:  $x\%$  (60% [P18], 75% [P19] or 85% [P20])  $\times$  (220 – your age).

After starting the programme, the resistance level will be controlled automatically constantly during exercise. You may also input target values for the function values time, distance and calories.

#### **P21 to P24 - User defined programmes (U1 to U4)**

In this mode you can exercise according to your own programmed profiles. After selecting a corresponding mode, you can preset each segment load in the profile (totally 10 segments) with the "UP" and "DOWN" buttons ("ENTER" for confirmation). The profile will be stored with the corresponding user-ID (U1...U4) and may be recalled later for another workout. After you have finished your input, press "ST/SP" to begin exercise. You may also input target values for different function values.

#### *Hint:*

- *The pulse values registered and displayed with this computer are merely indicative in nature and can deviate considerably from the actual pulse.*

# Using a chest-belt

The computer is equipped with a built-in receiver for the wireless receipt of the pulse frequency with a compatible SKANDIKA chest-belt (frequency: 5.3 kHz, not included in the scope of delivery). Recommended range: 100 - 150 cm.

If you want to transmit your pulse frequency via chest-belt, follow these instructions:

1. Make sure, the correct battery is inserted into the chest belt (button cell of type CR-2032, 3 volts).  
Battery replacement: Open the battery case cover on the back side of your transmitter (you may use a coin for easier handling) and replace the old battery by a fitting new one. The „+“ symbol should point upwards. Make sure, the sealing ring ist correctly placed. Close the case cover properly.
2. The transmitter needs to be attached to the belt with the 2 buttons. Moisten the electrodes slightly before wearing the belt. This increases the contact sensibility and enables a better wireless signal transmitting.
3. Put on the chest belt and make sure the electrodes are touching your breast directly. You may adjust the belt to a comfortable length. However, the electrodes must have contact with your skin in order to measure your heart rate correctly.
4. The 5.3 kHz-chest-belt communicates with the computer via the built-in receiver. It may last up to one minute before a correct display can be seen.

## **Safety hint:**

Do not use the chest-belt if you have a pacemaker, unless you double-checked with the producer of the pacemaker and with your doctor!

## *Battery-Disposal:*

Batteries should not be considered as regular garbage. As consumer you are obliged to return finished batteries. The finished batteries can be returned to a collection base at your residential area or at places, where batteries can be bought.

# Using the apps „Fitness Data“ & „KinoMap“

The computer can be connected via **Bluetooth®** to communicate with the apps Fitness Data or KinoMap.

## 1) Fitness Data

Activate **Bluetooth®** on your mobile device, download the app and install on your iOS- or Android device. You may use the app search function of your mobile device or the following QR code to find the app:

Android System

iOS System

After successful installation you can start the app on your mobile device. First enter all personal data and settings (language, nickname, sex, unit of measure, age, height, weight). You may also upload a photo by clicking on the picture. Click on the green button on the bottom part of the screen („register“ or „sign up“). You may now enter target values on the next app screen or navigate to further pages:

Tap on bike icon to automatically search for available **Bluetooth®** computers of exercise devices. If the computer of your bike is not found, check whether the **Bluetooth®** function is not activated properly or if the computer is in sleep mode. If the computer is in sleep mode, start to pedal on your bike to wake it up. Reinitiate the automatic search.

If the **Bluetooth®** computer has been found, a prompt will appear. We recommend to search for the computer manually if the automatic search is not successful. You may search for an available computer to connect with the app on the page „device management“.

If the computer has been connected successfully once, the following connections can be established automatically in most cases. Click on bike icon and make sure, the computer is not in sleep mode. Automatic connection mode will be started if no operation follows for approx. 5 seconds.

After the app has been successfully connected, it will show the current training data:

Switch to selection page of exercise devices

Multiple function value display for time, speed, distance, odometer, calories, pulse etc. If „SCAN“ is activated, all available function values will be shown in sequence for approx. 6 seconds each

Status of the **Bluetooth®**-connection:  
Green = Connection o.k.  
White = No connection to the computer device

Shows the available function values of your training session

Tap on „History“, to recall the records of the last 30 days. Here you can see different statistics of your training, e.g. the achievement degree. You may also select single days to get more information on single training sessions.

On the page „Music Player“ you may select local stored music files for playback:

Playback random music files

Mute

Repeat one or all music files


Slide to the left to show music file list

## 2) KinoMap

The free KinoMap app turns your device into a powerful and fun fitness product! Choose from the geotagged video on Kinomap and go! Your position is shown in the video and as icon on the map. Download the Kinomap Fitness app in the app or play store by searching for 'Kinomap Fitness' or follow one of the links below.

Apple App store: <https://itunes.apple.com/nl/app/kinomap-fitness/id611177969?mt=8>

Google Play Store: <https://play.google.com/store/apps/details?id=com.kinomap.fitness>

Install and open the app and register resp. log in. Activate Bluetooth®. Now add the exercise device as follows:

**1.** Tap the menu icon.

**2.** Tap 'Equipment' in the sidebar.

**3. („+“)** Tap the '+' button in the top right corner of the 'Equipment' screen.

**4. („Exercise bike“)** Tap 'Exercise bike' in the list of device types.


**5. („Skandika“)** Tap the 'Skandika' brand logo in the selection grid.


**6. Tap on the device** Tap on the device name in the 'Discovered' list.

Tap on OK and return to main menu.

You can now choose a course for your exercise (depending on your type of KinoMap account maybe only some free courses). Tap on „Start“ and begin!

# Exploded diagram

## Parts list

No.	Description	Q'ty	No.	Description	Q'ty
1	Main frame	1	37	Cover R	1
2	Front post	1	38	Sensor cable	1
3	Seat tube	1	39	Transformer wheel end (L &R)	2
4	Sliding bracket	1	40	Front stabilizer cover	1
5	Handlebars	1	41	Rear stabilizer cover R	1
6	Front stabilizer	1	42	Rear stabilizer cover L	1
7	Rear stabilizer	1	43	Front clamp cover	1
8	Seat (OPTIONAL)	1	44	Bearing	2
9	Computer	1	45	Spring	2
10	Hand pulse	2	46	Bushing	1
11	Foam grip	2	47	Cable	1
12	Cap	2	48	Upper tension cable	1
13	Rear clamp cover	1	49	Lower tension cable	2
14	Front clamp cover	1	50	Hand pulse wire	2
15	Motor	1	51	Screw M3*20	24
16	T-knob	1	52	Screw M4.2*15	14
17	End cap	2	53	Screw M4.2*20	1
18	Brake disc	1	54	Hexagon nut	1
19	Knob	1	55	Lower cable wire	5
20	Rear cover	1	56	Nylon nut	14
21	Cup	1	57	Flat washer	8
22	Release knob	1	58	Special bolt	4
23L	Pedal L	1	59	Pan bolt M8x55	4
23R	Pedal R	1	60	Adjustable bolt	2
24	Crank cover	2	61	Hexagon nut	1
25L	Crank L	1	62	Bracket	2
25R	Crank R	1	63	Stop nut	2
26	Chain cover	2	64	Stop nut	4
27	Cover L	1	65	Hexagon bolt	4
28	Rear stabilizer end cap (L & R)	2	66	Flat washer	4
29	Rear stabilizer cover	1	67	Spring shim	4
30	Flywheel	1	68	Nylon nut	4
31	Idle wheel with double bearing	1	69	Pan bolt	2
32	Idle wheel bracket	1	70	Flat washer	6
33	Spring	1	71	Pan bolt M8x16	6
34	Pulley	1	72	Bottle holder	1
35	Pulley axle	1	73	Adaptor	1
36	Belt	1			

## Cleaning and care

- Before cleaning the appliance, make sure that the appliance is switched off.
- Clean the device with a slightly damp cloth only. Never use brushes, aggressive cleaning agents, petrol, thinners or alcohol. Rub the device dry with a clean dry cloth.
- Wait until the device is completely dry before using it again.

## Guarantee conditions

For our devices we provide a warranty as defined below.

1. In accordance with the following conditions (numbers 2-5) we repair defect or damage to the device free of charge, if the cause is a manufacturing defect. Therefore, these defects / damages need to be reported to us without delay after appearance and within the warranty period of 24 months after delivery to the end user. The warranty does not cover parts, which easily break (e. g. glass or plastic). The warranty does not cover slight deviations of the product, which are insignificant for usability and value of the device and damage caused by chemical or electrochemical effects and damages caused by penetration of water or generally force majeure damage.
2. The warranty achievement is the replacement or repair of defective parts, depending on our decision. The cost of material and labor will be borne by us. Repairs at customer site cannot be demanded. The proof of purchase along with the date of purchase and / or delivery is required. Replaced parts become our property.
3. The warranty is void if repairs or adjustments are made, which are not authorized by us or if our devices are equipped with additional parts or accessories that are not adapted to our devices. Furthermore, the warranty is void if the device is damaged or destroyed by force majeure or due to environmental influences and in case of improper handling / maintenance (e.g. due to non-observance of the instruction manual) or mechanical damages. The customer service may authorize you to replace or repair defective parts after telephone consultation. In this case, the warranty is not void.
4. Warranty services do not extend the warranty period nor do they initiate a new warranty period.
5. Further demands, especially claims for damages which occurred outside the device, are excluded as long as a liability is not obligatory legal.
6. Our warranty terms - which cover the requirements and scope of our warranty conditions - do not affect the contractual warranty obligations of the seller.
7. Parts of wear and tear are not included in the warranty.
8. The warranty is void if not used properly or if used in gyms, rehabilitation centers and hotels. Even if most of our units are suitable for a professional use, this requires a separate agreement.

### Environmental protection

At the end of its life cycle, this product must not be disposed of with household waste but must be taken to a collection unit for the recycling of electric and electronic equipment. The symbol on the product, the instructions for use or the packaging express mention of this. The basic materials can be recycled as specified on the labelling.

When recycling the materials and finding other utilisation for used equipment, you are making a significant contribution towards protecting our environment. Ask at your council about the respective local disposal sites.

For service, accessories and spare parts, please contact: [info@skandika.com](mailto:info@skandika.com)

Service centre: MAX Trader GmbH, Wilhelm-Beckmann-Straße 19, 45307 Essen, Germany

In accordance with our policy of continual product improvement, we reserve the right to make technical and visual changes without notice.

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